Existing feature modification or enhancements:

- → platform/ui/src/components/Header/Header.tsx
 - Logo Updated:

The OHIF logo was **replaced** with the new AIH logo.

```
src={aihLogo}
    alt="Logo"
    style={{
        width: '45px',
        height: '45px',
        borderRadius: '50%',
        marginLeft: '10px',
        marginRight: '5px',
        objectFit: 'contain',
}}
```

• Local storage is cleared when the components mount:

```
useEffect(()=>{
   localStorage.setItem('items', '');
   window.dispatchEvent(new Event('itemsData'));
}, []);
```

• Replaced Icons.ByName and Icons.ChevronClosed with Icon:

```
<Icon
    name="chevron-left"
    className="w-8 text-white"
/>
```

• Patient data anonymisation:

```
showPatientInfo = PatientInfoVisibility.VISIBLE_COLLAPSED,showPatientInfo = PatientInfoVisibility.DISABLED,
```

- → platform/ui/tailwind.config.js
 - Font size updates:

```
- xxs: '0.625rem', // 10px
- xs: '0.6875rem', // 11px
- sm: '0.75rem', // 12px
- base: '0.8125rem', // 13px
- lg: '0.875rem', // 14px
- xl: '1rem', // 16px
- // 2xl and above will be updated in an upcoming version
+ xs: '0.65rem',
+ sm: '0.75rem',
+ base: '0.8125rem',
+ lg: '0.9275rem',
+ xl: '1.25rem',
```

• Colour Changed:

From blue to green

New Features:

→ Login Page :

platform/app/src/components/LoginSignup/LoginSignup.tsx

- Users can enter their email and password to log in.
- On successful login, authentication tokens are stored in the browser, and the user is redirected to the home page.
- Shows an error message if login credentials are invalid.

→ platform/ui/src/api/apiClient.js

Delete Study:

• Share Study:

```
+ async updateAccess(studyUid, email, role) {
+ ///console.log(studyUid);
+ const requestBody = {
+ email: email,
+ role: role,
+ };
```

```
+ ///console.log('sending share api',
apiEndpoints.updateAccess(studyUid));
+ const res = await this.makeRequest(
+ 'POST',
+ apiEndpoints.updateAccess(studyUid),
+ requestBody,
+ {},
+ true
+ );
+ ///console.log('RESS', res);
+ return res;
+ }
```

• Labeling Studies:

• Fetching study access list:

```
+ async getAccessList(studyUid) {
+ ///console.log('sending share list api',
apiEndpoints.getAccessList(studyUid));
+ return await this.makeRequest('GET',
apiEndpoints.getAccessList(studyUid), undefined, {}, true);
+ }
```

• Handle Models:

```
+ async handleModel(model, series = 'all', frames = 'all',
studyUid, setToastMessage) {
+ ...
+ }
```

```
+ async handleMammoModel(studyUid, setToastMessage) {
    ...
}
```

```
+async handleGBCModel(studyUid, setToastMessage) {
...
}
```

```
+ async handleXRayModel(studyUid, setToastMessage) {
...
}
```

→ Side Panel:

platform/ui/src/components/SidePanelModified/SidePanelModified.tsx

- Dropdown: Lets users switch between content for "Mammo," "GBC," and "X-Ray."
- Dynamically displays content based on the selected option.
- → Bi-Rads Scoring Feature:

extension/measurement-tracking/src/panels/MammoPanel/ActionButtons.tsx

→ GB-Rad Scoring Feature:

extensions/measurement-tracking/src/panels/GbcPanel/ActionButtons.tsx